



Compact and Affordable 8 Mph Scooter



PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

(Version 1: April 2017)

Features and Specification

Introduction

The compact Colt Sport Scooter features a sleek, sporty design and delivers high-performance operation. The standard features include a handbrake, front and rear suspension, five piece disassembly, a wraparound delta tiller, high-visibility automatic brake lights and an LED kerb light. These features combine to make the Colt Sport a superb scooter at an exceptional price.



Delta tiller control



Front and rear LED lights



Front and rear suspension

Available on



Available Colours

The Colt Sport comes in two optional colours:



Red



Blue

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
25 Stone (158.7 kg)	8 Mph (12.9 km/h)	5.7 cm (2.25 inch's)	133.35 cm (52.5 inch's)	119.4 cm (47 inch's)	65.7 cm (25.9 inch's)	Front: 25cm Rear: 25cm	10.7 Stone 68 k.g.	29.48 k.g. Front Section	Hi-Back Seat 20k.g	x2 40 Ah Battery

Warranty

1. Life time limited warranty on frame
2. Two year limited warranty on drive train, motor and brake
3. One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.

Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST —