

Jazzy air



The Jazzy Air's advanced suspension provides unique stability in raised and lower positions. Users are able to drive safely at 3.5 mph fully elevated enabling enhanced **#socialmobility**



...A new era in powered mobility.



Jazzy POWER CHAIRS



Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —



@PRIDEMOBILITYUK



Crash Tested

Tel: 01869 324600 | Email: sales@pride-mobility.co.uk | Web: www.pridemobility.co.uk

(Version 5: June 2018)

Introduction

Introducing the new Jazzy Air® a new era in power mobility.
The advanced stability and performance of the chassis allows the Jazzy Air to be safely elevated to 10", whilst driving at 3.5 mph.

The Jazzy Air® features Active-Trac® suspension, which provides unique stability in both the raised and lowered position. With a tight turning radius of 57.15 cm (22.5"), unique Mid-Wheel 6° Drive Design and the addition of bright LED markers combines to make the Air a smoother more enjoyable journey.

#socialmobility

Jazzy Air | Seating, Colour Options & Tech Spec

Seating Options

Choice of 16" x 16", 18" x 18" and 20" x 20" contoured high-back reclining sport seat

Colour Options



VR2 drive control with LED marker switch & Buddy button for 10" seat lift



10" Seat lift, which can be operated at 3.5 Mph











Integrated footplate and footplate design



Bright LED marker lights

Technical Specifications

Drive Wheels (cm)  25.04 cm solid tyres only	Caster Wheels (cm)  Front and Rear: 15.25 cm	Max Speed (mph)  4 mph (6.4 km/h)	Ground Clearance (cm)  2.54 cm At motor bracket	Turning Radius (cm)  57.15 cm	Overall Length (cm)  107.2 cm with foot rigging	Overall Width (cm)  64.77 cm with armrests	Seating Options (inches)  16" W - 16" D 18" W - 18" D 20" W - 20" D	Batteries  x 2 volt 30-35 AH Deep Cycle	Available Electronics  70A PG VR2	Weight Capacity (stone)  21 Stone 6 lbs (137 k.g)
Base Weight (kg)  82 kg	Battery Weight (kg)  11.1 kg	Warranty Lifetime limited warranty on frame. 1-year limited warranty on electronics. 18-month warranty on drive control system. 2-year limited warranty on drivetrain.				Notes 1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. 2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. 3. Battery weight may vary based on manufacturer.				