




*Exceptional
Turning Radius
and Dual Motor
Technology*

The Jazzy® Zero Turn features a 3wheel turning radius combined in a 4 wheel package. Featuring Pride's **NEW ZERO TURN TECHNOLOGY** (patent-pending), 21 Ah batteries, USB charging facilities, front and rear LED lights and much more.



@Pridemobilityuk



Version 1.5 | November 2017



32 Wedgwood Road, Bicester, Oxon, OX26 4UL | Tel: 01869 324600 | Fax: 01869 323070 |
Web: pride-mobility.co.uk | Email: sales@pride-mobility.co.uk



Jazzy Zero Turn | Standard Features



USB Charging Port



38.2" Turning Radius



Front and Rear Suspension



Front and Rear LED Lights



Under-seat Storage



21Ah Batteries


Jazzy Zero Turn | Colour Options



Fire Opal



Onyx Black



Pearl White



Pewter





Sapphire Blue



Green Moonstone

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
										
21 Stone 6lbs (136 kg)	4 Mph (6.4 km/h)	8.38 cm ⁴ (3.3 inch's)	97.2 cm (38.2 inch's)	116.2 cm (45.75 inch's)	50.8 cm (20 inch's)	Front: 17.78 cm Rear: 22.86 cm	8.5 Stone 54 k.g.	21.7 k.g. rear section	Folded Seat 12 k.g	x2 21 Ah deep cycle

Warranty

1. Life time warranty on frame
2. Two-year limited warranty on drive train, including: transaxle, motor and brake.
3. One-year limited warranty on all electrical parts
4. 12 month warranty for batteries

Notes

1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
3. Battery weight may vary based on manufacturer.
4. Measured at motor
5. Information correct at time of printing
6. Images used for illustration purposes only.